

To Wear Elaré is to live beyond time
youth. elegant. desired.

 ElaréTM
YOUTHFUL BEAUTY

USER MANUAL



Please read this user manual carefully before beginning and keep it properly for future reference
Before using this product for the first time, please perform a light sensitivity test.
Please refer to the instructions for details.

THE SPIRITUAL MEANING OF ELARÉ

Born from the essence of light itself, Elaré is a symbol of radiance, renewal, and timeless beauty. It represents the quiet glow that lives within every woman—the light that cannot be measured by age, only by the way she shines.

Like the first rays of dawn after a long night, Elaré embodies hope, confidence, and transformation. It is a reminder that beauty is not something to chase, but something to awaken. With every moment of self-care, the Elaré LED Red Light Therapy Mask becomes more than a skincare ritual—it becomes a sacred pause, a return to yourself, a celebration of the woman you are and the woman you continue to become.

Elaré is for those who believe that true beauty deepens with self-love, that elegance is timeless, and that youth is not an age, but an energy. It is the glow of feeling desired, the confidence of feeling seen, and the power of embracing your most luminous self.


To wear Elaré is to live beyond time—forever youthful, effortlessly elegant, endlessly desired.

Catalogue

Product Description	1
Operating Mode	2
Function	3
Product Controller	4
How To Use	5
Specifications	5
Installation	6
Safety Cautions	7
Warning	8

Product Description

The wearable LED light therapy silicone facial mask is a household wearable and rechargeable photon rejuvenation beauty device. The surface structure of this product is composed of medical silicone material, which is convenient and easy to use. It is foldable and portable, and it can be fixed with the adjustable buckle. It has a total of 186 high-quality LEDs (face 122 LEDs, neck 64LEDs), a total of 744chips. Each LED emits 4 different wavelengths (Blue light 415nm, Yellow light 590nm, Red light 660nm, and Near-infrared light 850nm), and can be adjusted by the controller to emit 3 different energies. Based on the theory of phototherapy, the light wave stimulates collagen, elastin, fights acne causing bacteria, improve pigmentation and skin tone whilst soothing inflammation, calming and healing the skin.




Controller



Silicone safety goggles



Carrying Bag



Type-C cable




Adjustable buckle




User Manual

Modes



M1

Increases collagen, elastin and blood circulation




M2

Increases collagen, elastin, blood circulation, evens skin tone and helps reduce pigmentation



M3

Increases collagen, elastin and blood circulation. Calms redness and irritation.




M4

A combination of all LED modes except blue light



M5

Fights acne causing bacteria



M6

A combination of all LED modes

Function

6 Modes

- M1: Increases collagen, elastin and blood circulation ■
- M2: Increases collagen, elastin, blood circulation, evens skin tone and helps reduce pigmentation ■ ■
- M3: Increases collagen, elastin and blood circulation. Calms redness and irritation ■ ■
- M4: A combination of all LED modes except blue light ■ ■ ■
- M5: Fights acne causing bacteria ■
- M6: A combination of all LED modes ■ ■ ■ ■ ■



Yellow Light

Balances skin texture, stimulates red blood cells, reduces redness, swelling, and fine lines




Red Light

Anti-aging, promoting collagen, tightening and increasing elasticity



Near-infrared Light

Repair skin, improve blood circulation, and increase skin moisture




Blue Light

Sterilizes acne producing bacteria, reduces inflammation and scars, and increases skin oxygen

Tips

Before opening the packaging for use, be sure to fully charge the battery of the mask. The charging time is about 3 hours. The upper cover is equipped with an adjustable buckle, allowing you to wear the mask comfortably. To install the mask, thread the strap through the entrance hole on the other side of the mask, wrap the strap around the head and firmly secure it.

Product Controller



How To Use


- | Power button: Press and hold for about 3 seconds, the light will turn on.
- | The brightness selection key: short press to adjust and switch energy output mode 50% 75% 100%.
- | Timing: The default timing for each mode is 10 minutes when the mode is selected.
- | Mode selection: Press the power button after the light is turned on. There are 6 modes you can choose.
- | Pulse selection: you can press the pulse function key in any mode.

Product Specifications


Product Name	Wearable LED Light Therapy silicone Facial Mask		
Charging Voltage	5V/2A	Battery capacity	5000mAh
Working Current	≤ 3.5A	Working rated power	≤ 15W
Face Piece	310*207*4.5mm	Neck Piece	350*260*4.5mm
LED Quantity	Face Piece 122 LEDs Neck Piece 64 LEDs	Wavelengths	415nm:590nm:660nm:850nm 1:1:1:1

Installation


1. Remove dustproof film









2. Put on the eye protector
(Silicone safety goggles)



3. Put on the adjustable
buckle



4. Connect power and controller



Safety Cautions

Light Sensitivity Test

NOTE: DO NOT USE THIS PRODUCT IF YOU ARE ONE OF THE FOLLOWING PEOPLE

1. People who are pregnant, preparing for pregnancy, or breastfeeding;
2. people with light sensitivity or taking any medication that may cause photosensitivity; people with hereditary diseases of the eye (such as lupus erythematosus, photosensitive eczema or albinism, etc.)

Step 1: Place the device on your arm and turn on the switch to irradiate for 10 minutes. If symptoms such as redness or itching appear on the skin, it indicates that you may have photosensitive skin and it is recommended that you stop using it, or consult a dermatologist before using it again.

Step 2: Before use, please clean your face and keep it dry (without makeup) to prevent photosensitive ingredients that may be contained in cosmetics from affecting the efficacy of skin care.

Step 3: After use, the skin may appear red, which is a normal reaction of the skin, and will return to normal within 1-2 hours. If you feel any other discomfort during the process, please stop using it immediately.

Warnings

1. Please strictly follow the correct method of using this product and do not use the equipment outside of the instructions, otherwise, it may cause machine malfunction or reaction after.
2. Do not use this product to improve any other condition except for indications.
3. Pregnant women, breastfeeding women, and children should consult doctor
4. If the user suffers from physical, sensory, or is undergoing dermatology disease improvement, including the onset period (such as lupus erythematosus, photosensitive eczema, or albinism), it is not recommended to use this product. If it must be used, it must be under the supervision of a specially assigned person, or under the guidance of a responsible person (such as a dermatologist).
5. If you suffer from any photosensitivity disorder (photosensitivity), do not use this product. If you suffer from photosensitivity disorders, you may cause serious skin reactions.
6. Do not use on the same day as skincare products containing photosensitive ingredients.
7. Pay attention to protecting your eyes during use, wear eye protectors correctly.
8. If there is any discomfort during the use of this product, please stop using it immediately.
9. This product is not waterproof. Do not immerse the equipment in water or rinse it with water. If you think it has come into contact with water, do not use it.
10. Do not use it for infants or children, and store it out of reach of children.
11. If any parts are damaged, please do not use this product.
12. This product is only for daily use. Do not disassemble, modify, or repair it without authorization.
13. It is recommended to use this product 3-4 times a week, with a single-use time of 10-20 minutes.
14. 12 months guarantee from date of purchase.



Manufacturer:

A&A Mangement
9485-1136 Québec inc.

Address:

6179 av. Durocher, Outremont (Québec), H2V 3Y7 Canada

Email:

Support@elareskincare.com